

## TRI-HARD Three Part Cold Mounting Compound



MetPrep's New TRI-HARD Cold Mounting compound is a three part acrylic system comprising a powder and two liquids.

For use where maximum hardness and excellent edge retention is required. Cure time 15 minutes

11 10 41	TRI-HARD Kit - Regular 1.5Kg TRI-HARD Powder	Mixing Instructions
	1.0 litre TRI-HARD Liquid 2 0.5 litre TRI-HARD Liquid 1	Add Liquid 1 to Liquid 2 and stir thoroughly
11 10 49	TRI-HARD Kit - Large 6.0 Kg TRI-HARD Powder 4.0 litre TRI-HARD Liquid 2 2.0 litre TRI-HARD Liquid 1	Add the liquid mixture to the powder and stir until homogeneous keeping bubbles to a minimum.

Mixing Ratio: Three parts Powder - Two parts Liquid 2 - One part Liquid 1

MetPrep Technical Data Sheet 0004 - 2015. For further information please contact sales@metprep.co.uk

MetPrep Ltd - Curriers Close, Charter Avenue, Coventry, CV4 8AW - 02476 421222