



## QUICK-SET 2 Ultra Fast Setting Acrylic Mounting Compound



MetPrep's Ultra Fast powder / liquid acrylic mounting system offers a hard filled mounting compound that exhibits minimal shrinkage with low odour,

A quality ultra fast setting option for use when specimen mounting speed is critical

11 10 35	QUICK-SET 2 powder - Hard Blue	1 Kg
11 10 36	QUICK-SET 2 liquid	0.5 Litre

Curing time 5 minutes

Mixing Instructions: Two Parts Powder - One Part Liquid

MetPrep Technical Bulletin 0006 - 2015. For further information please contact [sales@metprep.co.uk](mailto:sales@metprep.co.uk)

MetPrep Ltd - Curriers Close, Charter Avenue, Coventry, CV4 8AW - 02476 421222